

Product datasheet for **TP725451**

Pit1 (POU1F1) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human PIT-1(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	G61-R291
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human PIT-1 is produced by E.coli expression system and the target gene encoding G61-R291 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	5449
UniProt ID:	P28069
Synonyms:	Pituitary-specific positive transcription factor 1; Growth hormone factor 1 ; GHF1; POU1F1
Protein Families:	Druggable Genome, Transcription Factors



[View online »](#)