

## Product datasheet for **TP725461**

### RTEL1 Human Recombinant Protein

#### Product data:

|                                       |  |
|---------------------------------------|--|
| Product Type:                         | Recombinant Proteins   |
| Description:                          | Recombinant Human Regulator of telomere elongation helicase 1/RTEL1(C-His)   |
| Species:                              | Human  |
| Expression cDNA Clone or AA Sequence: | K800-T1000   |
| Tag:                                  | C-His  |
| Buffer:                               | Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4  |
| Note:                                 | MADE TO ORDER   Recombinant Human Regulator of telomere elongation helicase 1/RTEL1 is produced by E.coli expression system and the target gene encoding K800-T1000 is expressed with a His tag at the C-terminus. |
| Stability:                            | 12 months from date of despatch  |
| Locus ID:                             | 51750  |
| UniProt ID:                           | <a href="#">Q9NZ71</a>   |



[View online »](#)