

## Product datasheet for **TP725469**

### **FUBP1 Human Recombinant Protein**

#### **Product data:**

<b>Product Type:</b>	Recombinant Proteins
<b>Description:</b>	Recombinant Human Far upstream element-binding protein 1/FUBP1(C-His)
<b>Species:</b>	Human
<b>Expression cDNA Clone or AA Sequence:</b>	D341-S612
<b>Tag:</b>	C-His
<b>Buffer:</b>	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
<b>Note:</b>	MADE TO ORDER   Recombinant Human Far upstream element-binding protein 1/FUBP1 is produced by E.coli expression system and the target gene encoding D341-S612 is expressed with a His tag at the C-terminus.
<b>Stability:</b>	12 months from date of despatch
<b>Locus ID:</b>	8880
<b>UniProt ID:</b>	<a href="#">Q96AE4</a>



[View online »](#)