

## Product datasheet for **TP725478**

### Histone acetyltransferase MYST3 (KAT6A) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Histone acetyltransferase KAT6A(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	E800-S1000
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER   Recombinant Human Histone acetyltransferase KAT6A is produced by E.coli expression system and the target gene encoding E800-S1000 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	7994
UniProt ID:	<a href="#">Q92794</a>



[View online »](#)