

## Product datasheet for **TP725482**

### **KDM5C Human Recombinant Protein**

#### **Product data:**

<b>Product Type:</b>	Recombinant Proteins
<b>Description:</b>	Recombinant Human Lysine-specific demethylase 5C/KDM5C(C-His)
<b>Species:</b>	Human
<b>Expression cDNA Clone or AA Sequence:</b>	D1300-L1560
<b>Tag:</b>	C-His
<b>Buffer:</b>	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
<b>Note:</b>	MADE TO ORDER   Recombinant Human Lysine-specific demethylase 5C/KDM5C is produced by E.coli expression system and the target gene encoding D1300-L1560 is expressed with a His tag at the C-terminus.
<b>Stability:</b>	12 months from date of despatch
<b>Locus ID:</b>	8242
<b>UniProt ID:</b>	<a href="#">P41229</a>



[View online »](#)