

Product datasheet for **TP725493**

SPEN Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Msx2-interacting protein(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	V2-G271
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Msx2-interacting protein is produced by E.coli expression system and the target gene encoding V2-G271 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	23013
UniProt ID:	Q96T58
Synonyms:	SMART/HDAC1-associated repressor protein; SPEN homolog
Protein Families:	Druggable Genome



[View online »](#)