

## Product datasheet for **TP725497**

### MUS81 Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Crossover junction endonuclease MUS81(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	R202-Q420
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER   Recombinant Human Crossover junction endonuclease MUS81 is produced by E.coli expression system and the target gene encoding R202-Q420 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	80198
UniProt ID:	<a href="#">Q96NY9</a>



[View online »](#)