

Product datasheet for **TP725498**

Cleavage stimulation factor 2 (CSTF2) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Cleavage stimulation factor subunit 2/CSTF2(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	A2-A248
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Cleavage stimulation factor subunit 2/CSTF2 is produced by E.coli expression system and the target gene encoding A2-A248 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	1478
UniProt ID:	P33240



[View online »](#)