

## Product datasheet for **TP725499**

### **HOXA5 Human Recombinant Protein**

#### **Product data:**

<b>Product Type:</b>	Recombinant Proteins
<b>Description:</b>	Recombinant Human Homeobox protein Hox-A5(C-His)
<b>Species:</b>	Human
<b>Expression cDNA Clone or AA Sequence:</b>	F64-P270
<b>Tag:</b>	C-His
<b>Buffer:</b>	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
<b>Note:</b>	MADE TO ORDER   Recombinant Human Homeobox protein Hox-A5 is produced by E.coli expression system and the target gene encoding F64-P270 is expressed with a His tag at the C-terminus.
<b>Stability:</b>	12 months from date of despatch
<b>Locus ID:</b>	3202
<b>UniProt ID:</b>	<u><a href="#">P20719</a></u>



[View online »](#)