

Product datasheet for **TP725512**

HOXB3 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Homeobox protein Hox-B3(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	G177-Q396
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Homeobox protein Hox-B3 is produced by E.coli expression system and the target gene encoding G177-Q396 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	3213
UniProt ID:	P14651
Synonyms:	Homeobox protein Hox-2.7; Homeobox protein Hox-2G



[View online »](#)