

## Product datasheet for **TP725513**

### NEI3 (NEIL3) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Endonuclease 8-like 3/NEIL3(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	C257-Q490
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER   Recombinant Human Endonuclease 8-like 3/NEIL3 is produced by E.coli expression system and the target gene encoding C257-Q490 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	55247
UniProt ID:	<a href="#">Q8TAT5</a>
Synonyms:	DNA glycosylase FPG2; DNA glycosylase/AP lyase Neil3; Endonuclease VIII-like 3; Nei-like protein 3
Protein Families:	Druggable Genome
Protein Pathways:	Base excision repair



[View online »](#)