

## Product datasheet for **TP725514**

### ELF5 Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human ETS-related transcription factor Elf-5(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	L21-V245
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER   Recombinant Human ETS-related transcription factor Elf-5 is produced by E.coli expression system and the target gene encoding L21-V245 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	2001
UniProt ID:	<a href="#">Q9UKW6</a>
Synonyms:	E74-like factor 5; Epithelium-restricted ESE-1-related Ets factor; Epithelium-specific Ets transcription factor 2; ESE2
Protein Families:	Transcription Factors



[View online »](#)