

Product datasheet for **TP725536**

SNRPN Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human snRNP-N(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	T2-M215
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human snRNP-N is produced by E.coli expression system and the target gene encoding T2-M215 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	6638
UniProt ID:	P63162
Synonyms:	Small nuclear ribonucleoprotein-associated protein N; Sm protein D; Sm-D; Sm protein N; Sm-N; SmN; Tissue-specific-splicing protein
Protein Families:	Stem cell - Pluripotency



[View online »](#)