

Product datasheet for **TP725541**

POLR3F Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human DNA-directed RNA polymerase III subunit RPC6(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	A2-V260
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human DNA-directed RNA polymerase III subunit RPC6 is produced by E.coli expression system and the target gene encoding A2-V260 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	10621
UniProt ID:	Q9H1D9



[View online »](#)