

Product datasheet for **TP725546**

MAZ Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Myc-associated zinc finger protein/MAZI(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	A166-V420
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Myc-associated zinc finger protein/MAZI is produced by E.coli expression system and the target gene encoding A166-V420 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	4150
UniProt ID:	P56270



[View online »](#)