

Product datasheet for **TP725547**

HEY1 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human HEY1 Protein(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	S38-V240
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human HEY1 Protein is produced by E.coli expression system and the target gene encoding S38-V240 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	23462
UniProt ID:	Q9Y5J3
Synonyms:	Hairy/enhancer-of-split related with YRPW motif protein 1; Cardiovascular helix-loop-helix factor 2; CHF-2;Class B basic helix-loop-helix protein 31; bHLHb31; HES-related repressor protein 1; Hairy and enhancer of split-related protein 1; HESR-1;Hairy-related transcription factor 1; HRT-1; hHRT1
Protein Families:	Druggable Genome, Transcription Factors



[View online »](#)