

Product datasheet for **TP725551**

Histone H1.1 (HIST1H1A) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Histone H1.1(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	A11-K215
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Histone H1.1 is produced by E.coli expression system and the target gene encoding A11-K215 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	3024
UniProt ID:	Q02539



[View online »](#)