

Product datasheet for **TP725553**

MEF2D Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Myocyte-specific enhancer factor 2D/MEF2D(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	N293-K521
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Myocyte-specific enhancer factor 2D/MEF2D is produced by E.coli expression system and the target gene encoding N293-K521 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	4209
UniProt ID:	Q14814



[View online »](#)