

## Product datasheet for **TP725554**

### MEF2A Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Myocyte-specific enhancer factor 2A/MEF2A(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	G2-V216
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER   Recombinant Human Myocyte-specific enhancer factor 2A/MEF2A is produced by E.coli expression system and the target gene encoding G2-V216 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	4205
UniProt ID:	<a href="#">Q02078</a>



[View online »](#)