

Product datasheet for **TP725560**

TACC3 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human TACC3 Protein(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	E689-I838
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human TACC3 Protein is produced by E.coli expression system and the target gene encoding E689-I838 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	10460
UniProt ID:	Q9Y6A5
Synonyms:	Transforming acidic coiled-coil-containing protein 3; ERIC-1



[View online »](#)