

Product datasheet for **TP725565**

FOXO3 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Forkhead box protein O3/FOXO3(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	A373-G673
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Forkhead box protein O3/FOXO3 is produced by E.coli expression system and the target gene encoding A373-G673 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	2309
UniProt ID:	O43524



[View online »](#)