

Product datasheet for **TP725574**

PPM1B Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Protein phosphatase 1B(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	G2-G378
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Protein phosphatase 1B is produced by E.coli expression system and the target gene encoding G2-G378 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	5495
UniProt ID:	O75688
Synonyms:	Protein phosphatase 2C isoform beta; PP2C-beta
Protein Families:	Druggable Genome, Phosphatase, Stem cell - Pluripotency
Protein Pathways:	MAPK signaling pathway



[View online »](#)