

## Product datasheet for **TP725576**

### **CXCL17 Human Recombinant Protein**

#### **Product data:**

<b>Product Type:</b>	Recombinant Proteins
<b>Description:</b>	Recombinant Human 6-Cys CXCL17(C-His)
<b>Species:</b>	Human
<b>Expression cDNA Clone or AA Sequence:</b>	L24-L119
<b>Tag:</b>	C-His
<b>Buffer:</b>	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
<b>Note:</b>	MADE TO ORDER   Recombinant Human 6-Cys CXCL17 is produced by E.coli expression system and the target gene encoding L24-L119 is expressed with a His tag at the C-terminus.
<b>Storage:</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
<b>Stability:</b>	12 months from date of despatch
<b>Locus ID:</b>	284340
<b>UniProt ID:</b>	<a href="#">Q6UXB2</a>
<b>Synonyms:</b>	C-X-C motif chemokine 17; Dendritic cell and monocyte chemokine-like protein; DMC; VEGF coregulated chemokine 1; 4-Cys CXCL17



[View online »](#)