

Product datasheet for **TP725580**

HDGF Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human HDGF(P201K)(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	S2-L240(P201K)
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human HDGF is produced by E.coli expression system and the target gene encoding S2-L240(P201K) is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	3068
UniProt ID:	P51858



[View online »](#)