

## Product datasheet for **TP725584**

### LSP1 Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Lymphocyte-specific protein 1/LSP1(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	A2-P339
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER   Recombinant Human Lymphocyte-specific protein 1/LSP1 is produced by E.coli expression system and the target gene encoding A2-P339 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	4046
UniProt ID:	<u><a href="#">P33241</a></u>



[View online »](#)