

## Product datasheet for **TP725593**

### KCNIP4 Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human KCHIP4(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	N2-I250
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER   Recombinant Human KCHIP4 is produced by E.coli expression system and the target gene encoding N2-I250 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	80333
UniProt ID:	<a href="#">Q6PIL6</a>
Synonyms:	Kv channel-interacting protein 4; A-type potassium channel modulatory protein 4; Calsenilin-like protein; Potassium channel-interacting protein 4
Protein Families:	Druggable Genome, Ion Channels: Other



[View online »](#)