

## Product datasheet for **TP725595**

### **BIN2 Human Recombinant Protein**

#### **Product data:**

|  |  |
|--|--|
| <b>Product Type:</b>                         | Recombinant Proteins   |
| <b>Description:</b>                          | Recombinant Human Bridging integrator 2/BIN2(C-His)  |
| <b>Species:</b>                              | Human  |
| <b>Expression cDNA Clone or AA Sequence:</b> | A2-N244  |
| <b>Tag:</b>                                  | C-His  |
| <b>Buffer:</b>                               | Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4  |
| <b>Note:</b>                                 | MADE TO ORDER   Recombinant Human Bridging integrator 2/BIN2 is produced by E.coli expression system and the target gene encoding A2-N244 is expressed with a His tag at the C-terminus.                     |
| <b>Storage:</b>                              | Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles. |
| <b>Stability:</b>                            | 12 months from date of despatch  |
| <b>Locus ID:</b>                             | 51411  |
| <b>UniProt ID:</b>                           | <a href="#">Q9UBW5</a>   |
| <b>Synonyms:</b>                             | Breast cancer-associated protein 1   |



[View online »](#)