

## Product datasheet for **TP725598**

### Myosin Phosphatase (PPP1R12A) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Protein phosphatase 1 regulatory subunit 12A/PPP1R12A(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	A707-E981
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER   Recombinant Human Protein phosphatase 1 regulatory subunit 12A/PPP1R12A is produced by E.coli expression system and the target gene encoding A707-E981 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	4659
UniProt ID:	<a href="#">O14974</a>



[View online »](#)