

Product datasheet for **TP725600**

CoCoA (CALCOCO1) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Calphoglin(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	L341-E691
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Calcium-binding and coiled-coil domain-containing protein 1 is produced by E.coli expression system and the target gene encoding L341-E691 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	57658
UniProt ID:	Q9P1Z2
Synonyms:	Calphoglin; Coiled-coil coactivator protein; Calcium-binding and coiled-coil domain-containing protein 1; Sarcoma antigen NY-SAR-3; CALCOCO1
Protein Families:	Transcription Factors



[View online »](#)