

Product datasheet for **TP725613**

AMPK beta 1 (PRKAB1) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human AMPK subunit beta-1(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	G2-I270
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human AMPK subunit beta-1 is produced by E.coli expression system and the target gene encoding G2-I270 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	5564
UniProt ID:	Q9Y478



[View online »](#)