

Product datasheet for **TP725617**

ANKRD54 Human Recombinant Protein

Product data:

| | |
|--|---|
| Product Type: | Recombinant Proteins |
| Description: | Recombinant Human Ankyrin repeat domain-containing protein 54/ANKRD54(C-His) |
| Species: | Human |
| Expression cDNA Clone or AA Sequence: | R201-R300 |
| Tag: | C-His |
| Buffer: | Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4 |
| Note: | MADE TO ORDER Recombinant Human Ankyrin repeat domain-containing protein 54/ANKRD54 is produced by E.coli expression system and the target gene encoding R201-R300 is expressed with a His tag at the C-terminus. |
| Storage: | Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles. |
| Stability: | 12 months from date of despatch |
| Locus ID: | 129138 |
| UniProt ID: | Q6NXT1 |
| Synonyms: | Lyn-interacting ankyrin repeat protein |



[View online »](#)