

## Product datasheet for **TP725621**

### LRP1B Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human LRP-1B(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	L21-E244
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER   Recombinant Human LRP-1B is produced by E.coli expression system and the target gene encoding L21-E244 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	53353
UniProt ID:	<a href="#">Q9NZR2</a>
Synonyms:	Low-density lipoprotein receptor-related protein 1B;Low-density lipoprotein receptor-related protein-deleted in tumor; LRP-DIT
Protein Families:	Druggable Genome, Transmembrane



[View online »](#)