

Product datasheet for **TP725623**

BNIP3L Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human BCL2/BNIP3L Protein(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	S2-E184
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human BCL2/BNIP3L Protein is produced by E.coli expression system and the target gene encoding S2-E184 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	665
UniProt ID:	O60238
Synonyms:	adenovirus E1B 19 kDa protein-interacting protein 3-like; Adenovirus E1B19K-binding protein B5; BCL2/adenovirus E1B 19 kDa protein-interacting protein 3A; NIP3-like protein X; NIP3L
Protein Families:	Druggable Genome, Transmembrane



[View online »](#)