

Product datasheet for **TP725632**

RPA34 (RPA2) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human RP-A p32(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	L141-E270
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human RP-A p32 is produced by E.coli expression system and the target gene encoding L141-E270 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	6118
UniProt ID:	<u>P15927</u>



[View online »](#)