

## Product datasheet for **TP725633**

### FGF12 Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human FGF-12(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	E67-T243
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER   Recombinant Human FGF-12 is produced by E.coli expression system and the target gene encoding E67-T243 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	2257
UniProt ID:	<a href="#">P61328</a>
Synonyms:	Fibroblast growth factor 12; Fibroblast growth factor homologous factor 1;FHF-1; Myocyte-activating factor
Protein Families:	Secreted Protein
Protein Pathways:	MAPK signaling pathway, Melanoma, Pathways in cancer, Regulation of actin cytoskeleton



[View online »](#)