

Product datasheet for **TP725634**

Liver Carboxylesterase 1 (CES1) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Liver carboxylesterase 1/CES1(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	S385-E563
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Liver carboxylesterase 1/CES1 is produced by E.coli expression system and the target gene encoding S385-E563 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	1066
UniProt ID:	P23141



[View online »](#)