

Product datasheet for **TP725641**

Hsp40 (DNAJB1) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Dnaj homolog subfamily B member 1/DNAJB1(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	S151-I340
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Dnaj homolog subfamily B member 1/DNAJB1 is produced by E.coli expression system and the target gene encoding S151-I340 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	3337
UniProt ID:	P25685



[View online »](#)