

Product datasheet for **TP725643**

KPC2 (UBAC1) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human UBA domain-containing protein 1(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	I64-T405
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human UBA domain-containing protein 1 is produced by E.coli expression system and the target gene encoding I64-T405 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	10422
UniProt ID:	Q9BSL1



[View online »](#)