

## Product datasheet for **TP725644**

### PRKRA Human Recombinant Protein

#### Product data:

|                                       |                                                                                                                                                                                                                                                                              |
|---------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type:                         | Recombinant Proteins                                                                                                                                                                                                                                                         |
| Description:                          | Recombinant Human PRKRA(C-His)                                                                                                                                                                                                                                               |
| Species:                              | Human                                                                                                                                                                                                                                                                        |
| Expression cDNA Clone or AA Sequence: | S2-H200                                                                                                                                                                                                                                                                      |
| Tag:                                  | C-His                                                                                                                                                                                                                                                                        |
| Buffer:                               | Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4                                                                                                                                                                                          |
| Note:                                 | MADE TO ORDER   Recombinant Human PRKRA is produced by E.coli expression system and the target gene encoding S2-H200 is expressed with a His tag at the C-terminus.                                                                                                          |
| Storage:                              | Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.                                                                 |
| Stability:                            | 12 months from date of despatch                                                                                                                                                                                                                                              |
| Locus ID:                             | 8575                                                                                                                                                                                                                                                                         |
| UniProt ID:                           | <a href="#">O75569</a>                                                                                                                                                                                                                                                       |
| Synonyms:                             | Interferon-inducible double-stranded RNA-dependent protein kinase activator A; PKR-associated protein X; PKR-associating protein X; Protein activator of the interferon-induced protein kinase; Protein kinase, interferon-inducible double-stranded RNA-dependent activator |
| Protein Families:                     | Druggable Genome                                                                                                                                                                                                                                                             |



[View online »](#)