

## Product datasheet for **TP725645**

### **TTC11 (FIS1) Human Recombinant Protein**

#### **Product data:**

<b>Product Type:</b>	Recombinant Proteins
<b>Description:</b>	Recombinant Human Mitochondrial fission 1 protein/FIS1(C-His)
<b>Species:</b>	Human
<b>Expression cDNA Clone or AA Sequence:</b>	E2-G122
<b>Tag:</b>	C-His
<b>Buffer:</b>	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
<b>Note:</b>	MADE TO ORDER   Recombinant Human Mitochondrial fission 1 protein/FIS1 is produced by E.coli expression system and the target gene encoding E2-G122 is expressed with a His tag at the C-terminus.
<b>Stability:</b>	12 months from date of despatch
<b>Locus ID:</b>	51024
<b>UniProt ID:</b>	<a href="#">Q9Y3D6</a>



[View online »](#)