

## Product datasheet for **TP725649**

### TRBP (TARBP2) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human RISC-loading complex subunit TARBP2(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	S2-L279
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER   Recombinant Human RISC-loading complex subunit TARBP2 is produced by E.coli expression system and the target gene encoding S2-L279 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	6895
UniProt ID:	<a href="#">Q15633</a>
Synonyms:	TAR RNA-binding protein 2; Trans-activation-responsive RNA-binding protein
Protein Families:	Stem cell - Pluripotency



[View online »](#)