

## Product datasheet for **TP725652**

### MSF (SEPT9) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Septin-9(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	G250-G568
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER   Recombinant Human Septin-9 is produced by E.coli expression system and the target gene encoding G250-G568 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	10801
UniProt ID:	<a href="#">Q9UHD8</a>
Synonyms:	MLL septin-like fusion protein MSF-A; MLL septin-like fusion protein; Ovarian/Breast septin; Ov/Br septin; Septin D1
Protein Families:	Druggable Genome



[View online »](#)