

Product datasheet for **TP725667**

S100 Calcium Binding Protein A13 (S100A13) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Protein S100-A13(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	A2-K98
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Protein S100-A13 is produced by E.coli expression system and the target gene encoding A2-K98 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	6284
UniProt ID:	Q99584



[View online »](#)