

Product datasheet for **TP725670**

PDZK1 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human NHERF-3(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	T2-M519
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human NHERF-3 is produced by E.coli expression system and the target gene encoding T2-M519 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	5174
UniProt ID:	Q5T2W1
Synonyms:	Na(+)/H(+) exchange regulatory cofactor NHE-RF3; CFTR-associated protein of 70 kDa; Na(+)/H(+) exchanger regulatory factor 3; Na/Pi cotransporter C-terminal-associated protein 1 (NaPi-Cap1); PDZ domain-containing protein 1; Sodium-hydrogen exchanger regulatory factor 3



[View online »](#)