

## Product datasheet for **TP725673**

### RBM5 Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human RNA-binding protein 5/RBM5(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	G2-C204
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER   Recombinant Human RNA-binding protein 5/RBM5 is produced by E.coli expression system and the target gene encoding G2-C204 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	10181
UniProt ID:	<a href="#">P52756</a>



[View online »](#)