

Product datasheet for **TP725677**

Acyltransferase like 1 (LPCAT2) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human LPCAT-2(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	E340-D544
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human LPCAT-2 is produced by E.coli expression system and the target gene encoding E340-D544 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	54947
UniProt ID:	Q7L5N7
Synonyms:	Lysophosphatidylcholine acyltransferase 2; LPC acyltransferase 2; LysoPC acyltransferase 2
Protein Families:	Transmembrane



[View online »](#)