

Product datasheet for **TP725681**

N WASP (WASL) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Actin nucleation-promoting factor WASL(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	F31-N234
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Actin nucleation-promoting factor WASL is produced by E.coli expression system and the target gene encoding F31-N234 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	8976
UniProt ID:	<u>O00401</u>
Synonyms:	Neural Wiskott-Aldrich syndrome protein; N-WASP
Protein Families:	Druggable Genome
Protein Pathways:	Adherens junction, Chemokine signaling pathway, Fc gamma R-mediated phagocytosis, Pathogenic Escherichia coli infection, Regulation of actin cytoskeleton



[View online »](#)