

## Product datasheet for **TP725686**

### TGN46 (TGOLN2) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human TGOLN2 Protein(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	S25-E381
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER   Recombinant Human TGOLN2 Protein is produced by E.coli expression system and the target gene encoding S25-E381 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	10618
UniProt ID:	<a href="#">O43493</a>
Synonyms:	Trans-Golgi network integral membrane protein 2; Trans-Golgi network glycoprotein 46; Trans-Golgi network glycoprotein 48;Trans-Golgi network glycoprotein 51; Trans-Golgi network protein 2; TGN38 homolog; hTGN46; hTGN48; hTGN51
Protein Families:	Transmembrane



[View online »](#)