

## Product datasheet for **TP725687**

### VEGFB Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human VEGF-B(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	P22-A207
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER   Recombinant Human VEGF-B is produced by E.coli expression system and the target gene encoding P22-A207 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	7423
UniProt ID:	<a href="#">P49765</a>



[View online »](#)