

## Product datasheet for **TP725688**

### MAVS Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human MAVS (C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	G50-T300
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER   Recombinant Human MAVS is produced by E.coli expression system and the target gene encoding G50-T300 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	57506
UniProt ID:	<a href="#">Q7Z434</a>



[View online »](#)